



# Smarter software for better dentures.


Dentsply Sirona Digital Dentures combine our time-tested and trusted materials with AvaDent's unmatched design expertise and intelligent software. Continuously optimized by computer-aided engineering, these dentures are designed and produced through a sophisticated set of algorithms based on tens of thousands of real cases generating a truly Digital Denture.



Preview of Digital Denture design in relation to reference denture with transparency overlay.



Virtual simulation of Digital Denture robotic milling process.



InLab design software simulation of Digital Denture with puck in MCX5 holder.

## Avadent Design Software

### Powerful Software

Intelligent software capable of converting traditional inputs into an accurately fitting denture customized for appearance and optimal performance.

### Monitoring in Real Time

A dynamic 3D preview function allows the dental professional to review the digital setup and request any adjustments prior to fabrication, ensuring better-looking and better-fitting dentures for the patient.

### A Lasting Record

Should the patient desire a duplicate or replacement denture at any time, the dentist can order an identical denture to be produced from the patient's unique Digital Denture record. The Teeth Forever® Card contains the essential case information to reorder duplicate dentures.